HEMOSTATIX®

Hemostatix Thermal Scalpel Handle

Model 5050

For use with the Hemostatix Thermal Scalpel System Model 150, 3000 and 600 series controllers

The Hemostatix Thermal Scalpel Handle is a surgical cutting instrument designed to retain the precise, clean cutting characteristics of a traditional steel scalpel and to minimize blood loss by simultaneously sealing blood vessels as they are cut. The Hemostatix Thermal Scalpel System provides hemostasis with minimum tissue damage and virtually no muscle stimulation, using heat thermally conducted to the tissue from an elevated-temperature blade.

The Hemostatix Thermal Scalpel System consists of three components:

1. **Controller.** An electronic power supply which energizes the blade and provides various automatic calibration, sensing, and control functions. It has user controls with visual and audible indications of instrument status.
2. **Handle.** A resposable scalpel handle, connected to the Controller by a lightweight flexible electrical cable, for use with the disposable scalpel blades.
3. **Blades.** Different sizes of sterile disposable scalpel blades are available (similar in size, shape, and sharpness to conventional blades), incorporating heating and microcircuitry.

HOW SUPPLIED

The Hemostatix Thermal Scalpel Handle is supplied sterile if the pouch is unopened and undamaged. Prior to any subsequent reuse, the handle must be resterilized. (Refer to HANDLE CLEANING AND RESTERILIZATION section.)

BLADES

The model 5050 handle is designed for use with HEMOSTATIX Scalpel Blades. 57XX Series Blades are used with the Model 600 and 3000 Controllers. 56XX Series Blades are used with the Model 150 Controller.

Hemostatix Medical Technologies • 8400 Wolf Lake Dr. #109 • Bartlett, TN 38133 U.S.A.
901.261.0012 • Fax 901.261.0017
For information: 1.800.SCALPEL
www.hemostatix.com
1-800-SCALPEL

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EU Rep: QFI Ltd., 11 Burford Road
Stratford, London E152ST
United Kingdom
0473

Controller Front Panel

- 1. Change Instrument Indicator (Model 600 only)**
Lights (red) if the blade is not properly seated, or if blade is outside of the calibration limit, or if there is no blade in the handle.
- 1. Change Blade Indicator (Model 150 and 3000 Controller only)**
Lights or flashes (red) if the blade is not properly seated, or if blade is outside of the calibration limit. On the Model 150 Controller, the Change Blade Indicator also flashes if there is no blade in the handle.
- 2. Calibrating Indicator**
Lights or flashes (yellow) while the system is calibrating to a blade. Turns off when calibration is complete.
- 3. Ready Indicator (Model 3000 Controller Only)**
Lights (green) after calibration is completed and remains on as long as blade and handle are acceptable for use or until the blade or handle is unplugged or the Controller is turned off.
- 3. Ready to Heat Indicator (Model 150 Controller Only)**
Lights (yellow) after calibration is completed. Remains on until blade is activated by the handle On/Off switch.
- 4. Energized Indicator (Model 600 Only)**
Lights (amber) when blade is activated by the handle On/Off switch. Remains on as long as blade is activated.
- 4. Blade at Temp Indicator (Model 150 Controller Only)**
Lights (yellow) when blade is activated by the handle On/Off switch. Remains on as long as blade is activated.
- 5. Power Switch**
Turns on AC power.
- 6. Scalpel Connector**
Receptacle for connecting handle to Controller.
- 7. Volume Control**
Controls the audible tone volume continuously between the minimum and maximum levels.
- 8. Temperature Control**
Increases or decreases blade temperature.
- 9. Temperature Setting Display**
Displays blade temperature setting.

Handle Controls, Blade and Connector

10. On/Off Switch

Sliding the switch rearward mechanically latches the switch On and activates the blade. Red dot is visible when the switch is in the "ON" position. Sliding the switch forward, "OFF," deactivates the blade.

11. Thermal Coag Bar

When the On/Off switch is in the "ON" position, depressing the portion of the bar activates the Thermal Coag mode as long as the bar is depressed. Releasing the bar causes the temperature to revert the original temperature setting.

12. Temperature Control Switch

With the On/Off switch in the "OFF" position, depressing one of the Temperature Switches increases or decreases the temperature setting.

The Up Arrow increases the temperature setting and the Down Arrow decreases the temperature setting. The temperature setting continues to increase or decrease as long as the Temperature Control switch is depressed until the preset limit is reached. The temperature setting remains at the last setting when the switch is released.

(Note: When using the Hemostatix Scalpel with the Model 150 Controller, only the Up Arrow is active. Once the highest setting is reached, holding the switch depressed resets the temperature to the lowest setting and resumes the upward change sequence. Depressing the Down Arrow results in no change to the temperature setting.)

13. Blade

Insert connector end of blade into handle until it snaps into place and is firmly seated.

14. Handle Connector

Connector for attaching handle to Controller.

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


Model 3000

Model 600D

SYSTEM OPERATING INSTRUCTIONS-MODEL 3000

1. Attach the handle to the Controller by centering the arrow on the handle connector at the top of the front panel receptacle and plugging it into the receptacle.
2. Turn on the Controller by activating the POWER switch. If a blade is not in the handle, no lights will be present except for the setting indicator on the Model 3000 Controller. When a functional blade is inserted, the CALIBRATING light will turn on while the Controller is calibrating to the blade. The CHANGE BLADE light will activate whenever the continuous calibration function determines that the blade has become unacceptable for use. When the Controller has finished calibrating and determined the blade to be acceptable for use, the READY light will go on.
3. To insert new blade: insert the connector end of the blade into handle (sharp edge on side opposite On/Off switch) until it snaps into place and is firmly seated.
4. Adjust the SCAPEL SET POINT control to the initial setting desired by the surgeon. Initial skin incision should be made with blade cold (handle in off position). For other tissues, typical settings are between 180°C and 300°C, depending upon the surgeon's preference and the type of tissue and/or surgical procedure involved. The TEMPERATURE SETTING may also be changed by depressing either the UP or DOWN indicators on the Temperature Control Switch on the side of the handle when the handle On/Off switch is in the "OFF" position.
5. To heat the blade, slide On/Off switch rearward to the "ON" position to expose the red dot on the handle. An audible intermittent tone will sound while the blade is heating, followed by a steady tone once the blade is at the selected temperature.
6. To activate the Thermal Coag mode, depress the Thermal Coag bar down. The SCAPEL SET POINT will then change to 300°C. As long as the Thermal coag bar is depressed, the Controller will sound an intermittent tone. When the Thermal Coag switch is released, the Controller will return to the original set point and a steady tone.
7. To change the blade temperature, operate the SCAPEL SET POINT switch located on the face of the Controller. The blade temperature may also be changed by pressing the Temperature Control Up/Down switches on the handle when the handle On/Off switch is in the "OFF" position.
8. An adjustable-volume signal provides a steady tone when the blade is at the selected temperature. The audible tone is intermittent if the blade is not the selected temperature, as during warm-up, Thermal Coag, or undergoing excessive power dissipation.
9. If at any time the blade fails to have acceptable electrical characteristics as determined by the controller, the CHANGE BLADE light will go on, and a warning will sound.
10. To remove a blade from the handle, pull straight out.
CAUTION: If the blade is hot, use a sponge or a clamp to grasp the blade on the black portion of the blade.
11. If the blade becomes dislodged, cool the blade to room temperature before reinsertion.

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SYSTEM OPERATING INSTRUCTIONS-MODEL 150 CONTROLLER

1. Attach the handle to the Controller by aligning the arrow on the handle connector with the arrow on the front panel receptacle and plugging it into the receptacle.
2. Turn on the Controller by activating the POWER switch. If a blade is not in the handle, the CHANGE BLADE light will flash. When a functional blade is inserted, the CALIBRATING light will flash while the Controller is calibrating to the blade. The CHANGE BLADE light will reactivate whenever the continuous calibration function determines that the blade has become unacceptable for use. When the Controller has finished calibrating and determined the blade to be acceptable for use, the READY TO HEAT light will go on.
3. To insert new blade: insert the connector end of the blade into handle (sharp edge on side opposite On/Off switch) until it snaps into place and is firmly seated.
4. Adjust the TEMPERATURE SETTING control to the original setting desired by the surgeon. Initial skin incision should be made with blade cold (handle in Off position). For other tissues, typical settings are between 180°C and 260°C, depending upon the surgeon's preference and the type of tissue and/or surgical procedure involved. The TEMPERATURE SETTING may also be changed by depressing the UP indicator on the Temperature Control Switch on the side of the handle when the handle "OFF" switch is in the Off position. (Note: the temperature setting of the Model 150 Controller is upwardly adjustable only. Depressing the "Down" indicator on the handle will not affect the temperature setting.)
5. To heat the blade, slide the On/Off switch rearward to the "ON" position to expose the red dot on the handle. The READY TO HEAT light will go out and an audible intermittent tone will sound while the blade is heating; the tone will change to a steady tone and the BLADE AT TEMP light will go on when the blade reached the selected temperature.
6. To activate the Thermal Coag mode, depress the Thermal Coag bar down. The TEMPERATURE SETTING DISPLAY will then change to 270°C. As long as the Thermal Coag bar is depressed, the Controller will sound an intermittent tone. When the Thermal Coag switch is released, the Controller will return to the original set point and a steady tone.
7. To change the blade temperature, operate the TEMPERATURE SETTING switch on the face of the Controller. The blade temperature may also be changed by pressing the Up Temperature Control switch on the handle when the handle On/Off switch is in the "OFF" position.
8. An adjustable-volume signal provides a steady tone when the blade is at the selected temperature. The audible tone is intermittent if the blade is not at the selected temperature, as during warm-up, Thermal Coag, or when undergoing excessive power dissipation.
9. If at any time the blade fails to have acceptable electrical characteristics as determined by the controller, the CHANGE BLADE light will go on, and the READY TO HEAT light or the BLADE AT TEMP light will go off.
10. To remove a blade from the handle, pull straight out.
CAUTION: if the blade is hot, use a sponge or a clamp to grasp the blade on the black portion of the blade.
11. If the blade becomes dislodged, cool the blade to room temperature before reinsertion.

SYSTEM OPERATING INSTRUCTIONS - MODEL 600 CONTROLLER

See the Hemostatix Thermal Scalpel Operating Manual for detailed instructions on use of the Controller and System.

1. Attach the handle to the Controller by aligning the arrow on the handle connector with the arrow on the front panel receptacle and plugging it into the front panel receptacle.
2. Apply power to the instrument by placing the POWER switch in the up position. After completion of the power on self test, the CHANGE INSTRUMENT light will illuminate if a blade is not in the handle. When a functional blade is inserted, the CHANGE INSTRUMENT light will extinguish, and the CALIBRATING light will illuminate. The CHANGE INSTRUMENT light will reactivate and an audible alarm is heard whenever the automatic and continuous error checking function determines that the blade has become unacceptable for further use.
3. Insert new blade. Insert connector end of blade into handle (sharp edge on side opposite On/Off switch) until the blade is firmly seated.
4. After successful calibration is complete the CALIBRATING light will go out. An audible success tone also indicates completion of the calibration.
5. Adjust the TEMPERATURE SETTING to the initial setting desired by the surgeon. Initial skin incision should be made with blade cold (handle in off position). For other tissues, typical settings are between 180°C and 300°C depending upon the surgeon's preference and the type of tissue and/or surgical procedure involved.
6. To heat the blade, slide On/Off switch to the "ON" position to expose the red dot. The ENERGIZED light will turn on and the audible tone is heard.
7. To activate the Thermal Coag mode, depress the Thermal Coag bar down. The TEMPERATURE SETTING DISPLAY will then change to 300°C. As long as the Thermal Coag bar is depressed, the Controller will sound an intermittent tone. When the Thermal Coag switch is released, the Controller will return to the original set point and a steady tone.
8. To change to blade temperature operate the TEMP switches UP or DOWN. The blade temperature may also be changed by pressing the Temperature Control switches on the handle when the Handle On/Off switch is in the "OFF" position.
9. An adjustable-volume signal provides a steady tone when the blade is at the selected temperature. The audible signal is intermittent if the blade edge is not at the selected temperature (as during warm-up or cool down), Thermal Coag, or excessive dissipation.
10. If at any time the blade fails to have acceptable electrical characteristics as determined by the instrument, the CHANGE INSTRUMENT light will illuminate, the ENERGIZED light will extinguish, and an audible warning tone will be generated.
11. To remove a blade from the handle, pull straight out.
CAUTION: If the blade is hot, use a sponge or a clamp to grasp the blade on the back portion of the blade.
12. If the blade becomes dislodged, cool the blade to room temperature before reinsertion.

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Handle Cleaning and Resterilization

Cleaning

DO NOT IMMERSE THE HANDLE IN LIQUID OF ANY KIND. Do not allow any solution to penetrate to the interior of the handle. Do not use ultrasonic cleaning.

Manual cleaning is the recommended method of cleaning for the Model 5050 handle and cable. Since the device cannot be immersed, it is important to adhere to the following guidelines.

Cleaning should begin as soon after use as possible. Prepare an enzyme solution (e.g. Klenzyme™) and a detergent solution (e.g., Manu-Klenz™), as recommended by the manufacturers. An enzyme is not a cleaner. An enzyme is intended to break down blood and bodily fluids (i.e., protein) and thus facilitate cleaning. Detergents contain wetting and emulsifying agents that suspend soil and prevent the formation of insoluble compounds or scum on the device or on the surface of the cleaning solution.

Gently wipe the surface of the device with a sponge or towel that has been moistened with the enzyme solution. Rinse the blade receptacle of the device thoroughly with the enzyme solution and use a soft bristle brush to remove any gross soil, paying particular attention to crevices and other hard-to-clean areas. The use of pipe cleaners or cotton swabs should be avoided since they tend to fray, leaving particles inside the blade receptacle. Carefully examine the device for any adherent visible soil (e.g., blood, protein). Enzymes are themselves protein substances and must be thoroughly removed with a detergent. Repeat the steps of wiping the surface of the device and cleaning the blade receptacle with a soft bristle brush and the detergent solution. Rinsing with lukewarm tap water is necessary to remove all traces of detergents and extraneous debris.

Following cleaning, shake the device to remove excess water and wipe the handle and cable with a soft, dry cloth.

Sterilization

Before resterilization, inspect to ensure there are no breaks in the cable and no foreign material in the blade receptacle or the connector. Ensure that the handle On/Off switch is in the "ON" or "OFF" position (not in between). Loop the cord assembly into approximately seven inch diameter circles and place in appropriate packaging material.

The sterilization processes predominantly used in health care facilities are steam sterilization and ethylene oxide (EO) sterilization. Both methods are effective in sterilizing the Model 5050 handle and cable. EO sterilization is the preferred sterilization process for it will provide the longest useful life of the device. However, steam sterilization is cost-effective and poses no concerns about exposure of personnel or patients to toxic sterilant residuals. Steam sterilization may be used, although repeated exposures at temperatures of 132°C (270°F) in an "unwrapped" or "flash" condition may result in discoloration and/or deformation of the handle parts, cord, or connector, thus reducing the life of the product. A complete drying cycle should be performed to reduce/prevent eventual corrosion from the steam sterilization process. If steam autoclaving is to be used, the handles are to be "wrapped" and the recommended parameters are:

Gravity displacement steam sterilization

Temperature:	132°C (270°F)
Exposure time:	15 minutes
Drying time:	30 minutes

Temperature:	121°C (250°F)
Exposure time:	30 minutes
Drying time:	30 minutes

Prevacuum steam sterilization

Temperature:	132°C (270°F)
Exposure time:	4 minutes
Drying time:	30 minutes

The above guidelines and parameters have been properly tested to provide the best methods for cleaning and resterilization of your device. Should you have further questions, please contact us at 1-800-SCALPEL.

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PRECAUTIONS

Precautions for use of the Controller are listed on the unit. The following should be emphasized.

- DO NOT IMMERSE HANDLE IN LIQUIDS.
- DO NOT USE THIS PRODUCT IN THE PRESENCE OF FLAMMABLE GASES OR LIQUIDS.
- HANDLE MUST BE RSTERILIZED PRIOR TO REUSE. See section on Handle Cleaning and Resterilization.
- Plastic Skin Drapes
If a plastic skin drape is used, initial skin incision should be made cold (handle in the "OFF" position).
- Electrosurgical Units

DO NOT TOUCH HEMOSTATIX BLADE TO ELECTROSURGICAL (Bovie) TIP SINCE SIGNIFICANT DAMAGE TO THE HEMOSTATIX SCALPEL CONTROLLER WILL RESULT. Keep at least 1 cm between the Hemostatix blade and the electrosurgical tip.

- Avoiding Inadvertent Patient Burns

DO NOT REST THE HEMOSTATIX SCALPEL ON THE DRAPES OR ON THE PATIENT WHEN BLADE IS ENERGIZED. When energized, the blade is sufficiently hot that patient burns can result from contact.

For portions of procedures where the Hemostatix scalpel is not being used, it is recommended that the handle On/Off switch can be positioned "OFF."

- Handle Drying Following Autoclaving

Immediate use after autoclaving may result in a temporary malfunction due to moisture retained within the handle. Use of a drying cycle in the autoclave or other means of drying the handle prior to use can preclude such problems.

DANGER: POSSIBLE EXPLOSION HAZARD IF USED IN THE PRESENCE OF FLAMMABLE ANESTHETICS.

Hemostatix Thermal Scalpel handles are limited-use, reusable devices intended to be replaced after no more than 10 uses.


Contacts in the handle will wear and oxidize with repeated use. When handle contacts have thus degraded, the handle will no longer function properly and must be replaced with a new one.

USE OF A HANDLE WITH DEGRADED CONTACTS CAN RESULT IN BLADES RUNNING HOTTER OR COOLER THAN THE "TEMPERATURE SETTING" VALUE AND MAY ADVERSELY AFFECT SURGICAL PERFORMANCE.


Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or other licensed practitioner.


Hemostatix is a registered trademark of Hemostatix Medical Technologies.

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WARRANTY

Hemostatix Medical Technologies hereby warrants to the original purchaser that reasonable care has been used in the manufacture of the Hemostatix Scalpel Handle and that, when properly used, it will be free from defects in the material or workmanship for a period of sixty (60) days after the date of shipment from Hemostatix Medical Technologies. The sole and exclusive remedy with respect to any Hemostatix Scalpel Handle found within its warranty period not to meet these standards is that, after return to and examination by Hemostatix Medical Technologies, Hemostatix Medical Technologies will without charge at its option either repair or replace that portion of the Hemostatix Scalpel Handle found to be defective. This warranty shall not apply (a) if that portion of the Hemostatix Scalpel Handle has been repaired by anyone other than qualified service personnel approved by Hemostatix Medical Technologies or altered in any way which, in Hemostatix Medical Technologies's judgment, affects its usability or reliability; or (b) if the serial or lot number has been altered, effaced, or removed; or (c) if the fault has been caused by abnormal conditions of operation or misuse by permitting electrical contact with an active electrosurgical (Bovie) electrode. In such cases, Hemostatix Medical Technologies will notify the purchaser of Hemostatix Medical Technologies's determination, and repairs if requested will be billed at Hemostatix Medical Technologies's normal rates. Estimates will be submitted before any repair work is started, if requested.

Equipment being returned must obtain an RGA# from Hemostatix Medical Technologies before it is sent. Handles being returned must be cleaned, sterilized, and packaged in sterile packaging which verifies the sterility of the handle prior to return to Hemostatix Medical Technologies. Any handle not properly cleaned, sterilized, and packaged as described in this warranty will be disposed of and no warranty will be in effect.

THE FOREGOING EXPRESS WARRANTY, AS CONDITIONED AND LIMITED, IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESSED OR IMPLIED, BY OPERATION OF LAW OR OTHERWISE INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OR MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. HEMOSTATIX MEDICAL TECHNOLOGIES SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, EXPENSE, OR LIABILITY DIRECT OR INDIRECT WITH RESPECT TO THIS PRODUCT. HEMOSTATIX MEDICAL TECHNOLOGIES NEITHER ASSUMES NOT AUTHORIZES ANY OTHER PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THIS PRODUCT.

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This product is not to be used on patients with implanted devices such as pacemakers, defibrillators, or other electronic devices. The use of this product may interfere with the normal operation of these devices. Patients with implanted devices should consult their physician before using this product. The manufacturer is not responsible for any adverse effects caused by the use of this product on patients with implanted devices.

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